


Push-To-Set™ Vacuum Regulator Ordering Matrix


Instructions: Select the adapter, vacuum regulator, and fitting to form an 11-digit part number. To order individual adapters or fittings see Fittings/ Adapters section of the Ohio Medical Suction and Oxygen Therapy catalog (SOT32).

Example: 8701-1224-901 **8701** Diamond Adapter **1224** Standard Three Mode Continuous Regulator **901** Tubing Nipple


Adapters 87XX



Regulator 1XXX



Fittings 9XXX




















* These models contain no ferrous metal that is magnetically attractive.

Wall Outlet — Adapter Identification

For proper outlet and adapter identification use chart below. Note the type of connector required, not the name of the outlet manufacturer. Note particularly if the Chemetron (NCG) outlet (#4 and #5) has a round or rectangular keyway below the screw. All Chemetron (NCG) O2 keyways are round. Old (pre 1966) Chemetron (NCG) vacuum and air keyways are round but new ones are rectangle.

Wall Outlets

Adapters

1.		Diamond		Diamond
2.		DISS	 	DISS Hand-I-Twist (HIT) DISS Nut & Gland
3.		Schrader®		Schrader®
4.		Chemetron® (NCG) Rectangle Striker (current vacuum and air)		Chemetron® (NCG) Rectangle Striker (current vacuum)
5.		Chemetron® (NCG) Round Striker (old vacuum/air outlets pre 1966) Old & Current O2		Chemetron® (NCG) Round Striker (Old vacuum/air outlets pre 1966) Old & Current O2
6.		Puritan Bennett®		Puritan Bennett®
7.		O.E.S. Twist		O.E.S. Twist
8.		O.E.S. Med Star®		O.E.S. Med Star®

For more information, contact your Ohio Medical Sales Specialist or authorized agent.

© 2014 Ohio Medical Corporation. All Rights Reserved.

This document contains information that is proprietary and confidential to Ohio Medical Corporation. Use of this information is under license from Ohio Medical Corporation. Any use other than that authorized by Ohio Medical Corporation is prohibited.

Ohio Medical Corporation and the Ohio Medical Corporation logo are registered trademarks and Push-To-Set is a trademark of Ohio Medical Corporation.

Puritan Bennett is a registered trademark of Nellcor Puritan Bennett LLC.

Schrader is a registered trademark of Schrader-Bridgeport International, Inc.

Chemetron and Med Star are registered trademarks of Allied Healthcare Products, Inc.

