

instavac[®]C

MOBILE ASPIRATOR

DESCRIPTION AND APPLICATION

Instavac[®]C is a high vacuum & flow suction unit. It produces constant suction with vacuum & flow surpassing standard wall suction. Instavac[®]C is packaged in a compact, durable molded Polystyrene case. The heart of the Instavac[®]C is a fan cooled carbon rotary vane pump. Instavac[®]C is insulated to provide one of the quietest sounding suction pumps on the market today. It is an excellent choice for the O.R., surgicenters, and doctors offices.

FEATURES AND BENEFITS

- Constant suction providing 550 mmHg and over 38 lpm flow
- Quiet operation
- Easy to clean and process for quick turnarounds
- Fan cooled pump, mechanical safety overflow trap and exhaust port filter
- 3 year warranty

IDEAL FOR

- G.I. Labs, Operating Rooms and Surgicenters requiring high vacuum and flow rates with quiet performance.
- Surgical, scoping, and most specialty procedures.



PORTABLE SUCTION SYSTEM

TECHNICAL SPECIFICATIONS

Performance

Vacuum Range: 0 to 550mmHg
Free Air FLOW: 38 LPM minimum
Pump: fan cooled rotary carbon vane type

Collection Device

Canister: disposable plastic with mechanical shutoff
Capacity: 1200cc standard
Tubing: 14" with hydrophilic bacterial filter

Controls

Vacuum Regulator: rotary type on panel
Vacuum Gauge: Calibrated to 500 mm Hg

Electrical Requirements

AC: 120 volts, 60 Hz, 3A

Physical

Weight: 34 lbs
Dimensions: 28"H x 14"W x 16"D
46"H with back pole

Purchasing Specifications

Instavac®C Part Number: 757000

Shipping Information

Shipping weight:
48 lbs.

Package Size:
30" H x 24" W x 24" D

Packaging:
One per carton



WITH RESPECT TO ELECTRIC SHOCK,
FIRE AND MECHANICAL AND OTHER
HAZARDS ONLY IN ACCORDANCE WITH
UL 2601-1 AND CAN/CSA
C22.2 NO. 601-1 #9N56

AVAILABLE ACCESSORIES

PART #	DESCRIPTION
756390	Disposable Chest Drainage Unit, Basket
756150	Dual Vacuum Port Assembly
AI4120010	Disp. Collection Canister, 1200cc (case of 10)
AI4120025	Disp. Collection Canister, 1200cc (case of 25)
AI2025	Bottle Cap Bracket for Slide Mount
AI2026	Bottle Cap Bracket for Aspirator Mount
AI2300C	Cap & Float Assembly, Disposable (package of 25)
AI2311C	Turn Around Kit, Cap & Float, Tubing Filter (set of 10)
AI5576C	Filter, Disposable Exhaust (package of 25)
AI5785-14C	Filter, In-line suction w/14" Tubing (case of 50)
AI5785-18C	Filter, In-line suction w/18" Tubing (case of 50)



2017 Ohio Medical LLC. This document contains information that is proprietary and confidential to Ohio Medical LLC. Use of this information is under license from Ohio Medical LLC. Any use other than that authorized by Ohio Medical LLC is prohibited. Ohio Medical and the Ohio Medical logo are registered trademarks of Ohio Medical LLC.