

# MiniOX<sup>®</sup> 3000 Oxygen Analyzer

Ideal for continuous oxygen monitoring, the MiniOX 3000 is portable, easy to use, and designed for in-line oxygen monitoring in hospital, home ventilator, NICU, and anesthesia delivery equipment applications. Powered by one 9-volt battery, the MiniOX 3000 Oxygen Monitor provides up to 1500 hours of use. Each unit is shipped complete with sensor, coiled cable, tee adapter, and battery.

## Health-Care Applications:

- Respiratory therapy
- Neonatal care
- Anesthesia

## Features:

- High and low oxygen concentration alarms
- Automatic calibration
- Battery life indicators with alarms
- Sensor/cable disconnect alarm
- Alarm verification indicator
- Extra-large display window
- Integral tabletop stand

## Part Numbers

MiniOX 3000	814365
MiniOX 3000 (International)	711451


**Note:** MiniOX 3000 ships with coiled cable, tee adapter, retaining strap, sensor with deflector, and battery.

## Accessories

Mounting Bracket (pole)	474664
Mounting Bracket (wall)	10023945
Coiled Cable	472045
Tee Adapter	473021
Retaining Strap	634249
Deflector	803229



## Specifications

Range	0-100% O <sub>2</sub>
Resolution	0.1 O <sub>2</sub> Increments
Linearity	± 1%
Battery Indication	<p><b>Stage 1:</b> Warning that approximately six hours of monitoring life remain. “LO BAT” (Global version displays “”) appears on display; beep sounds every 30 seconds until battery is replaced or Stage 2 occurs.</p> <p><b>Stage 2:</b> Alarm that proper operation is not assured. Monitor shuts down; audible and visual alarms are activated.</p>
Power Supply	One 9-volt battery
Battery Lifetime	Approximately 1500 hours
Cable Length	10 feet
Product Dimensions	5.98" x 3.26" x 1.31" (15.18 x 8.28 x 3.32 cm)
Product Weight	9.2 oz.
Operating Temperature Range	32°-104° F (0°-40° C)
Storage Temperature Range	-40°-158°F (-40°-70° C)
Humidity	5-95% RH (non-condensing)
Alarms	Pulsating alarm / flashing LED
Shipping Dimensions	10.25" x 7.56" x 6" (26.03 x 19.20 x 15.24 cm)
Shipping Weight	1.7 lb. (0.77 kg.)

© 2017 Ohio Medical, LLC. All rights reserved.

This document contains information that is proprietary and confidential to Ohio Medical, LLC. Use of this information is under license from Ohio Medical, LLC. Any use other than that authorized by Ohio Medical, LLC is prohibited.

Ohio Medical, MiniOX, and the Ohio Medical logo are registered trademarks of Ohio Medical, LLC.

1111 Lakeside Drive | Gurnee, IL 60031-4099 | Phone: 1-866-549-6446 | Fax: 1-847-855-6300 | [www.ohiomedical.com](http://www.ohiomedical.com)

550626 (Rev. 1, 02/2017)