



## Materials

- All brass in high-pressure and relief-valve zones
- 20-micron, sintered-bronze inlet filter
- Teflon pressure regulating seat
- Viton and silicone o-rings
- FDA-recommended brass and viton seal washer

## Specifications

Flows: 0-25 LPM  
Connection: CGA 870  
Outlet Style: Tubing Nipple  
Height: 5.75"  
Diameter: 1.48"  
PTO: 2

## Ordering Information

### GR25-870VX-TPP

25LPM 870 Click Style Gas Regulator with 2 power take-offs and a tubing nipple

# 3 Port Gas Regulators

## Click Style

### Product Features

- Patented brass core in all areas that contact high-pressure oxygen
- Meets ASTM-G-175-03 promoted ignition tests
- Inlet pressure: 200-3,000 PSI
- Outlet pressure: 50 PSI
- CGA 870 connection
- Available with two PTO outlets (50 PSI output pressure)
- Five-year warranty
- 0-25 LPM, FLOW SETTINGS: .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25



Visit [www.amvex.com](http://www.amvex.com)  
for the Gas Regulator part  
configurator and much more!